# HVACO<sup>®</sup>*endo* gel

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**Stand der Information: 03.2021**

## PRODUCT INFORMATION

Post surgical tissue adhesion is a common consequence of surgical procedures, these post-surgical adhesions are due to the formation of areas of contact, made of fibrous tissue, between adjacent internal organs, which may induce different post operative symptoms such as pain, infertility, bowel obstruction, restriction of movement, and thus a second surgery for removing fibrous tissue would be needed. In accordance to excellent surgery practice and preventing formation of post surgical adhesions it is recommended to use **HYAcorp endo gel** which forms a physical barrier against the contact between adjacent tissues and remain on the site of application for a sufficiently long period to avoid the formation of adhesions during the critical healing process of injured tissues. The product is degraded completely within six months.

**HYAcorp endo gel** is a bio-degradable, sterile, transparent, high viscous gel obtained by crosslinking special designed hyaluronic acid, one of the main components of human connective tissue and of epithelial and mesothelial tissues and it perfectly adheres to the tissue surface and to the abdominal wall creating an anti adhesion barrier.

### Composition

1 ml **HYAcorp endo gel** contains:

Hyaluronic acid	2.0 mg
Hyaluronic acid cross-linked	20.0 mg
Sodium chloride	6.9 mg
Water for injection ad	1.0 ml

### Indications

**HYAcorp endo gel** is indicated for laparotomic or laparoscopic abdominal interventions, and for hysteroscopic interventions for reduction of post-surgical adhesion formation.

### Instructions for use

1. After having brought the product to room temperature, open the pouch and transfer the syringe into the sterile field, adopting the normal aseptic techniques used in the surgical theatres.
2. Remove the protective cap on the tip of the syringe and insert the cannula / applicator.
3. After aspiration of all irrigation fluid, apply the gel inside to completely cover all anatomical sites where adhesion prevention is desired, even if the thickness of the gel applied does not influence the efficacy of the product, it is recommended to apply a 1 – 2 mm thick layer of gel.
4. Do not irrigate after application of the product.

### Special instructions for pelvic-gynecological surgeries

1. Lift ovary away from pelvic sidewall and apply a single layer of gel to cover the ovarian fossa and posterior surface of the ovary.
2. Return ovary to normal anatomical position and apply a single layer of gel to cover the anterior portion of the ovary.

## After placement – all surgeries:

1. Do not apply irrigation fluid to the sites coated with **HYAcorp endo gel**.
2. Do not attempt to reposition gel once applied.
3. After completely covering the site once, do not apply additional gel to that site, if the site is not adequately covered by a single layer of gel, additional layers can be applied next to previously applied gel. Avoid applying additional layers of gel on top of one another.

## Cannula

Cannulas suitable for the surgical procedure shall be used. For laparotomic procedures 18 G needles are recommended and for laparoscopic procedures 230 mm cannulas. Cannulas can be attached via the integrated LuerLock connector.

## Contraindications

Known hypersensitivity to hyaluronic acid and Gram positive bacteria. As all other implant materials, the device must not be used in patients with active infection and/or contamination of the surgical site. Bowel resection and anastomoses: Due to a possible risk of impairment with regard to the healing process, the product must not be used in bowel resections and surgeries involving anastomotic sites.

## Storage

Do not freeze or expose to extreme heat. Store at room temperature (2 °C – 25 °C). Freezing of the product can result in severe bubbles formation. Do not use product in this case as the performance cannot be guaranteed.

## Warning and precautions

This product is to be used only under surgical aseptic conditions by a licensed medical practitioner. Once package is opened discard any remaining gel after use. Do not use if package is damaged.

With regard to resections involving anastomosis sites there are no clinical data available to prove the safety and performance of **HYAcorp endo gel**.

## Presentation

The product is designed to satisfy all needs of different surgical procedures utilized by surgeons and gynecologists:

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**BioScience GmbH**